

Anti-RGPD2/3/4/6/8 antibody (577-626 aa) (STJ95433)

STJ95433

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-RANBP2-like and GRIP domain-containing protein 8 and RANBP2-like and GRIP domain-containing protein 5/6 and RanBP2-like and GRIP domain-containing protein 4 and RanBP2-like and GRIP domain-containing protein 3 and OBSO
Applications	WB/ELISA/IHC
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-2000
Range	IHC-P 1:50-300 ELISA 2000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	653489 285190 727851
Gene Symbol	RGPD3 RGPD4 RGPD3_HUMAN RGPD4_HUMAN RGPD8_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human RGPD2/3/4/6/8 at the amino acid range 577-626
Immunogen Region	577-626 aa
Specificity	RGPD2/3/4/6/8 Polyclonal Antibody detects endogenous levels of RGPD2/3/4/6/8 protein.
Immunogen Sequence	



Western blot analysis of the lysates from Jurkat cells using RGPD2/3/4/6/8 antibody.

Western blot analysis of RGPD2/3/4/6/8 Antibody. The lane on the right is blocked with the RGPD2/3/4/6/8 peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081