

Anti-Rsk-1/2/3/4 antibody (191-240 aa) (STJ95548)

STJ95548

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ribosomal protein S6 kinase alpha-1 and Ribosomal protein S6 kinase alpha-3 and Ribosomal protein S6 kinase alpha-2 and Ribosomal protein S6 kinase alpha-6 (191-240 aa) is suitable for use in Western Blot and ELISA res
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID [6196](#)
[27330](#)
[6197](#)
[RPS6KA2](#)
[RPS6KA6](#)
[KS6A2_HUMAN](#)
[KS6A6_HUMAN](#)
[KS6A3_HUMAN](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human RSK1/2/3/4 at the amino acid range 191-240

Immunogen Region 191-240 aa

Specificity Rsk-1/2/3/4 Polyclonal Antibody detects endogenous levels of Rsk-1/2/3/4 protein.

Immunogen Sequence



Western blot analysis of lysates from Jurkat cells, using RSK1/2/3/4 Antibody. The lane on the right is blocked with the synthesized peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081