

Anti-Sox-8/9/17/18 antibody (116-165 aa) (STJ95735)

STJ95735

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Transcription factor SOX-8 and Transcription factor SOX-9 and Transcription factor SOX-17 and Transcription factor SOX-18 (116-165 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and E
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Monkey


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	IHC 1:100-1:300 ELISA 1:10000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID [6662](#)
[54345](#)
[30812](#)
[SOX9](#)
[SOX18](#)
[SOX9_HUMAN](#)
[SOX18_HUMAN](#)
[SOX8_HUMAN](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human SOX8/9/17/18 at the amino acid range 116-165
Immunogen Region 116-165 aa
Specificity Sox-8/9/17/18 Polyclonal Antibody detects endogenous levels of Sox-8/9/17/18 protein.
Immunogen Sequence



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081