

Anti-STK32C antibody (340-389 aa) (STJ95825)

STJ95825

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Serine/threonine-protein kinase 32C (340-389 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 IHC 1:100-1:300 ELISA 1:5000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	282974
Gene Symbol	STK32C
Uniprot ID	ST32C_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human STK32C at the amino acid range 340-389
Immunogen Region	340-389 aa
Specificity	STK32C Polyclonal Antibody detects endogenous levels of STK32C protein.
Immunogen Sequence	



Immunohistochemistry analysis of paraffin-embedded human breast cancer, using STK32C Antibody. The picture on the right is blocked with the STK32C peptide.



Western blot analysis of STK32C Antibody. The lane on the right is blocked with the STK32C peptide.