

Anti-TBX18 antibody (121-170 aa) (STJ95932)

STJ95932

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-T-box transcription factor TBX18 (121-170 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	9096
Gene Symbol	TBX18
Uniprot ID	TBX18_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human TBX18 at the amino acid range 121-170
Immunogen Region	121-170 aa
Specificity	TBX18 Polyclonal Antibody detects endogenous levels of TBX18 protein.
Immunogen Sequence	



Western blot analysis of the lysates from HT-29 cells using TBX18 antibody.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081