

## Anti-TCEAL3/5/6 antibody (141-190 aa) (STJ95939)

STJ95939

### GENERAL INFORMATION


|                          |                                                                                                                                                                                                                                                            |
|--------------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                                                                                                                         |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-Transcription elongation factor A protein-like 3 and Transcription elongation factor A protein-like 5 and Transcription elongation factor A protein-like 6 (141-190 aa) is suitable for use in Immunohistochemistry, Immun |
| <b>Applications</b>      | IHC/IF/ELISA                                                                                                                                                                                                                                               |
| <b>Host/Source</b>       | Rabbit                                                                                                                                                                                                                                                     |
| <b>Reactivity</b>        | Human/Rat/Mouse                                                                                                                                                                                                                                            |

### PRODUCT PROPERTIES

|                      |                                                                                                                       |
|----------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>     | Polyclonal                                                                                                            |
| <b>Clone ID</b>      |                                                                                                                       |
| <b>Concentration</b> | 1 mg/mL                                                                                                               |
| <b>Conjugation</b>   | Unconjugated                                                                                                          |
| <b>Purification</b>  | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution</b>      | IHC 1:100-1:300                                                                                                       |
| <b>Range</b>         | ELISA 1:40000<br>IF 1:50-200                                                                                          |
| <b>Formulation</b>   | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| <b>Isotype</b>       | IgG                                                                                                                   |
| <b>Storage</b>       | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |
| <b>Instruction</b>   |                                                                                                                       |

### TARGET INFORMATION

|                           |                                                                                                                                               |
|---------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Gene ID</b>            | <a href="#">158931</a><br><a href="#">340543</a><br><a href="#">85012</a>                                                                     |
| <b>Gene Symbol</b>        | <a href="#">TCEAL6</a><br><a href="#">TCEAL5</a><br><a href="#">TCAL6_HUMAN</a><br><a href="#">TCAL5_HUMAN</a><br><a href="#">TCAL3_HUMAN</a> |
| <b>Immunogen</b>          | The antiserum was produced against synthesized peptide derived from the human TCEAL3 at the amino acid range 141-190                          |
| <b>Immunogen Region</b>   | 141-190 aa                                                                                                                                    |
| <b>Specificity</b>        | TCEAL3/5/6 Polyclonal Antibody detects endogenous levels of TCEAL3/5/6 protein.                                                               |
| <b>Immunogen Sequence</b> |                                                                                                                                               |



Immunohistochemistry analysis of paraffin-embedded human breast carcinoma, using TCEAL3/5/6 Antibody. The picture on the right is blocked with the synthesized peptide.