

**Anti-TCP1 antibody (Internal) (STJ95948)**

STJ95948

**GENERAL INFORMATION**

|                          |                                                                                                                                                   |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Type</b>      | Primary antibodies                                                                                                                                |
| <b>Short Description</b> | Rabbit polyclonal antibody anti-T-complex protein 1 subunit alpha (Internal) is suitable for use in Western Blot and ELISA research applications. |
| <b>Applications</b>      | WB/ELISA                                                                                                                                          |
| <b>Host/Source</b>       | Rabbit                                                                                                                                            |
| <b>Reactivity</b>        | Human/Mouse/Rat                                                                                                                                   |

**PRODUCT PROPERTIES**

|                            |                                                                                                                       |
|----------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Clonality</b>           | Polyclonal                                                                                                            |
| <b>Clone ID</b>            |                                                                                                                       |
| <b>Concentration</b>       | 1 mg/mL                                                                                                               |
| <b>Conjugation</b>         | Unconjugated                                                                                                          |
| <b>Purification</b>        | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Dilution Range</b>      | WB 1:500-1:2000<br>ELISA 1:10000                                                                                      |
| <b>Formulation</b>         | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.                                               |
| <b>Isotype</b>             | IgG                                                                                                                   |
| <b>Storage Instruction</b> | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.                        |

**TARGET INFORMATION**

|                           |                                                                                   |
|---------------------------|-----------------------------------------------------------------------------------|
| <b>Gene ID</b>            | 6950                                                                              |
| <b>Gene Symbol</b>        | TCP1                                                                              |
| <b>Uniprot ID</b>         | TCPA_HUMAN                                                                        |
| <b>Immunogen</b>          | Synthesized peptide derived from the Internal region of human TCP-1 Alpha.        |
| <b>Immunogen Region</b>   | Internal                                                                          |
| <b>Specificity</b>        | TCP-1 Alpha Polyclonal Antibody detects endogenous levels of TCP-1 Alpha protein. |
| <b>Immunogen Sequence</b> |                                                                                   |