

Anti-Tubulin Alpha 1/3/4 antibody (238-287 aa) (STJ96143)

STJ96143

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Tubulin alpha-1A chain and Tubulin alpha-1B chain and Tubulin alpha-1C chain and OBSOLETE and Tubulin alpha-3E chain and Tubulin alpha-4A chain (238-287 aa) is suitable for use in Western Blot and ELISA research applic
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:10000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	84790 112714 7277 TUBA1C TUBA3E TBA1C_HUMAN TBA3E_HUMAN TBA4A_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the human TUBA1/3/4 at the amino acid range 238-287
Immunogen Region	238-287 aa
Specificity	Tubulin Alpha 1/3/4 Polyclonal Antibody detects endogenous levels of Tubulin Alpha 1/3/4 protein.
Immunogen Sequence	