

Anti-Zic1/2/3 antibody (316-365 aa) (STJ96310)

STJ96310

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Zinc finger protein ZIC 1 and Zinc finger protein ZIC 2 and Zinc finger protein ZIC 3 (316-365 aa) is suitable for use in Western Blot, Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	IHC 1:100-1:300 IF 1:200-1:1000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID [7547](#)
[7546](#)
[7545](#)
[ZIC3](#)
[ZIC2](#)
[ZIC3_HUMAN](#)
[ZIC2_HUMAN](#)
[ZIC1_HUMAN](#)

Immunogen The antiserum was produced against synthesized peptide derived from the human ZIC1/2/3 at the amino acid range 316-365


Immunogen Region 316-365 aa

Specificity Zic1/2/3 Polyclonal Antibody detects endogenous levels of Zic1/2/3 protein.

Immunogen Sequence



Western blot analysis of lysates from Jurkat and COLO cells, using ZIC1/2/3 Antibody. The lane on the right is blocked with the synthesized peptide.



Immunohistochemistry analysis of paraffin-embedded human brain tissue, using ZIC1/2/3 Antibody. The picture on the right is blocked with the synthesized peptide.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081