

Anti-CYP2A6/7/13 antibody (Internal) (STJ96458)

STJ96458

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Cytochrome P450 2A6 and Cytochrome P450 2A7 and Cytochrome P450 2A13 (Internal) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:40000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	1548 1549 1553 CYP2A6 CYP2A7 CP2A6_HUMAN CP2A7_HUMAN CP2AD_HUMAN
Immunogen	Synthesized peptide derived from the Internal region of human CYP2A6/7/13.
Immunogen Region	Internal
Specificity	CYP2A6/7/13 Polyclonal Antibody detects endogenous levels of CYP2A6/7/13 protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081