

Anti-Ub antibody (STJ96797)

STJ96797

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat


PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID [6233](#)
[7316](#)
[7314](#)
[RPS27A](#)
[UBC](#)
[RS27A_HUMAN](#)
[UBC_HUMAN](#)
[UBB_HUMAN](#)

Immunogen Synthesized peptide derived from the human Ub around the non-acetylation site of K63.
Immunogen Region
Specificity Ub Polyclonal Antibody detects endogenous levels of Ub protein only when acetylation at K63.
Immunogen Sequence



Western blot analysis of HepG2 cells using Ub Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
 St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081