

Anti-UBA1 antibody (91-140 aa) (STJ96858)

STJ96858

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ubiquitin-like modifier-activating enzyme 1 (91-140 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-1:2000 ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7317
Gene Symbol	UBA1
Uniprot ID	UBA1_HUMAN
Immunogen	The antiserum was produced against synthesized peptide derived from the N-terminal region of human UBA1 at the amino acid range 91-140
Immunogen Region	91-140 aa
Specificity	UBA1 Polyclonal Antibody detects endogenous levels of UBA1 protein.
Immunogen Sequence	



Western blot analysis of PC12 cells using UBA1 Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081