

Anti-Ub antibody (N-Term) (STJ96863)

STJ96863

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C (N-Term) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	


TARGET INFORMATION

Gene ID [7314](#)
[7311](#)
[6233](#)
[UBB](#)
[UBA52](#)
[UBB_HUMAN](#)
[RL40_HUMAN](#)
[RS27A_HUMAN](#)

Immunogen Synthesized peptide derived from the N-terminal region of human Ub.
Immunogen Region N-Term

Specificity Ub Polyclonal Antibody detects endogenous levels of Ub protein.

Immunogen Sequence



Western blot analysis of HepG2, HeLa cells using Ub Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081