

Anti-Di-Methyl-Histone H3-Lys27 antibody [3B12] (STJ97003)

STJ97003

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Di-Methyl-Histone H3.1 and Histone H3.2 and Histone H3.3-Lys27 is suitable for use in Western Blot and Immunoprecipitation research applications.
Applications	WB/IP
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3B12
Concentration	
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:1000-3000 IP 1:200
Formulation	Liquid in PBS pH7.4, 0.5% BSA, 0.02% Sodium Azide and 50% Glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	126961/333932/653604 3020/3021 H3C15.H3C14.H3C13 H3-3A.H3-3B H32_HUMAN H33_HUMAN H31_HUMAN <
Immunogen	Synthetic Peptide of Histone H3 (Di Methyl Lys27)
Immunogen Region	
Specificity	The antibody detects endogenous Histone H3 (di methyl K27) protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081