

Anti-Acetyl-Ub-Lys29 antibody (STJ97362)

STJ97362

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C-Lys29 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES


Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION


Gene ID [7311](#)
[6233](#)
[7316](#)
[UBA52](#)
[RPS27A](#)
[RL40_HUMAN](#)
[RS27A_HUMAN](#)
[UBC_HUMAN](#)

Immunogen Synthesized acetyl-peptide derived from the human Ub around the acetylation site of K29.

Immunogen Region Acetyl-Ub (K29) Polyclonal Antibody
Specificity Synthesized peptide derived from the human Ub around the acetylation site of K29.
Immunogen Sequence



Western blot analysis of HepG2 cells using Acetyl-Ub (K29) Polyclonal Antibody. Secondary antibody was diluted at 1:20000



Western blot analysis of HepG2 cells using Acetyl-Ub (K29) Polyclonal Antibody. Secondary antibody was diluted at 1:20000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.

St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081