

Anti-Acetyl-Ub-Lys48 antibody (STJ97364)

STJ97364

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Acetyl-Ubiquitin-60S ribosomal protein L40 and Ubiquitin-40S ribosomal protein S27a and Polyubiquitin-B and Polyubiquitin-C-Lys48 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	WB 1:500-1:2000
Range	ELISA 1:20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7314 6233 7311 UBB RPS27A UBB_HUMAN RS27A_HUMAN RL40_HUMAN
Immunogen	Synthesized acetyl-peptide derived from the human Ub around the acetylation site of K48.
Immunogen Region	
Specificity	Acetyl-Ub (K48) Polyclonal AntibodySynthesized peptide derived from the human Ub around the acetylation site of K48.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081