

Anti-Rubisco-Large Chain antibody [3G7] (STJ97527)

STJ97527

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Rubisco-Large Chain is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Plant

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3G7
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	WB 1:2000-5000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Protein
Immunogen Region	
Specificity	Rubisco (Large Chain) protein detects endogenous levels of Rubisco (Large Chain)
Immunogen Sequence	



Western blot analysis of Arabidopsis with Rubisco (Large Chain) Mouse mAb diluted at 1) 1:2, 000, 2) 1:5, 000

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081