

Anti-TF antibody [5F9] (STJ97600)

STJ97600

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Serotransferrin is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research applications.
Applications	IHC/IF/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	5F9
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using specific immunogen.
Dilution Range	IHC 1:100-200 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	7018
Gene Symbol	TF
Uniprot ID	TRFE_HUMAN
Immunogen	Protein
Immunogen Region	
Specificity	Transferrin protein detects endogenous levels of Transferrin
Immunogen Sequence	



Immunohistochemical analysis of paraffin-embedded rat-liver using antibody diluted at 1:50.



Immunohistochemical analysis of paraffin-embedded human-liver using antibody diluted at 1:50.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081