

Anti-FABP4 antibody (61-121 aa) [5H11A4E11] (STJ97815)

STJ97815

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Fatty acid-binding protein, adipocyte (61-121 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5H11A4E11
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	2167
Gene Symbol	FABP4
Uniprot ID	FABP4_HUMAN
Immunogen	Purified recombinant fragment of A-FABP (aa61-121) expressed in E. Coli.
Immunogen Region	61-121 aa
Specificity	A-FABP Monoclonal Antibody detects endogenous levels of A-FABP protein.
Immunogen Sequence	



Western blot analysis using A-FABP monoclonal antibody against truncated Trx-FABP4 recombinant protein (1).