

## Anti-ITGB2 antibody [10E12] (STJ97909)

STJ97909

### GENERAL INFORMATION


<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Integrin beta-2 is suitable for use in Immunofluorescence and ELISA research applications.
<b>Applications</b>	IF/ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	10E12
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	IF 1:200-1:1000 ELISA 1:10000
<b>Formulation</b>	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	3689
<b>Gene Symbol</b>	ITGB2
<b>Uniprot ID</b>	ITB2_HUMAN
<b>Immunogen</b>	Purified recombinant fragment of CD18 expressed in E. Coli.
<b>Immunogen Region</b>	
<b>Specificity</b>	CD18 Monoclonal Antibody detects endogenous levels of CD18 protein.
<b>Immunogen Sequence</b>	



Confocal immunofluorescence analysis of HL60 cells using CD18 monoclonal antibody (green). Blue: DRAQ5 fluorescent DNA dye.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081