

Anti-CDK1 antibody [5A6] (STJ97946)

STJ97946

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Cyclin-dependent kinase 1 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse/Rat/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5A6
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	983
Gene Symbol	CDK1
Uniprot ID	CDK1_HUMAN
Immunogen	Purified recombinant fragment of human Cdk1 expressed in E. Coli.
Immunogen Region	
Specificity	Cdk1 Monoclonal Antibody detects endogenous levels of Cdk1 protein.
Immunogen Sequence	



Western blot analysis using Cdk1 monoclonal antibody against HeLa (1), Jurkat (2), KS62 (3), A431 (4), MCF-7 (5), RAW264.7 (6), NIN-213 (7), PC-12 (8), and Cos7 (9) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081