

Anti-CCND3 antibody [6H4] (STJ97979)

STJ97979

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-G1/S-specific cyclin-D3 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	6H4
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	896
Gene Symbol	CCND3
Uniprot ID	CCND3_HUMAN
Immunogen	Purified recombinant fragment of human Cyclin D3 expressed in E. Coli.
Immunogen Region	
Specificity	Cyclin D3 Monoclonal Antibody detects endogenous levels of Cyclin D3 protein.
Immunogen Sequence	



Western blot analysis using Cyclin D3 monoclonal antibody against NIH/3T3 (1) and Jurkat (2) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081