

Anti-EPHA5 antibody (620-774 aa) [8B10B1/8B10F5] (STJ98027)

STJ98027

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ephrin type-A receptor 5 (620-774 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8B10B1/8B10F5
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2044
Gene Symbol	EPHA5
Uniprot ID	EPHA5_HUMAN
Immunogen	Purified recombinant fragment of EphA5 (aa620-774) expressed in E. Coli.
Immunogen Region	620-774 aa
Specificity	EphA5 Monoclonal Antibody detects endogenous levels of EphA5 protein.
Immunogen Sequence	



Western blot analysis using EphA5 monoclonal antibody against truncated EPHA5-His recombinant protein (1) and truncated EPHA5 (aa620-774) -hlgGfC transfected CHO-K1 cell lysate (2).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081