

Anti-ETV1 antibody (1-191 aa) [1C8B6] (STJ98041)

STJ98041

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-ETS translocation variant 1 (1-191 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1C8B6
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	2115
Gene Symbol	ETV1
Uniprot ID	ETV1_HUMAN
Immunogen	Purified recombinant fragment of ER81 (aa1-191) expressed in E. Coli.
Immunogen Region	1-191 aa
Specificity	ER81 Monoclonal Antibody detects endogenous levels of ER81 protein.
Immunogen Sequence	



Western blot analysis using ER81 monoclonal antibody against truncated Trx-ETV1 recombinant protein (1) and full-length ETV1 (aa1-477) -hlgGfC transfected CHO-K1 cell lysate (2).