

Anti-FOXD3 antibody [5G9] (STJ98084)

STJ98084

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Forkhead box protein D3 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human/Mouse/Monkey

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	5G9
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	27022
Gene Symbol	FOXD3
Uniprot ID	FOXD3_HUMAN
Immunogen	Purified recombinant fragment of human FoxD3 expressed in E. Coli.
Immunogen Region	
Specificity	FoxD3 Monoclonal Antibody detects endogenous levels of FoxD3 protein.
Immunogen Sequence	



Western blot analysis using FoxD3 monoclonal antibody against NTERA-1 (1), HUVE-12 (2), HEK293 (3), HeLa (4), Jurkat (5), NIH/3T3 (6), K562 (7), RAW264.7 (8) and COS7 (9) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081