

Anti-GATA5 antibody [1B11] (STJ98102)

STJ98102

GENERAL INFORMATION



Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Transcription factor GATA-5 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1B11
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	140628
Gene Symbol	GATA5
Uniprot ID	GATA5_HUMAN
Immunogen	Purified recombinant fragment of human GATA-5 expressed in E. Coli.
Immunogen Region	
Specificity	GATA-5 Monoclonal Antibody detects endogenous levels of GATA-5 protein.
Immunogen Sequence	



Western blot analysis using GATA-5 monoclonal antibody against HEK293 (1) and GATA5-hlgGfC transfected HEK293 (2) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081