

Anti-IL1A antibody [1F3B3/3A6F11] (STJ98170)

STJ98170

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Interleukin-1 alpha is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1F3B3/3A6F11
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG2a
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3552
Gene Symbol	IL1A
Uniprot ID	IL1A_HUMAN
Immunogen	Purified recombinant fragment of human IL-1a expressed in E. Coli.
Immunogen Region	
Specificity	IL-1 Alpha Monoclonal Antibody detects endogenous levels of IL-1 Alpha protein.
Immunogen Sequence	



Western blot analysis using IL-1 Alpha monoclonal antibody against truncated IL-1 Alpha recombinant protein.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081