

Anti-LIPA antibody [9G7F12/7G6D7] (STJ98212)

STJ98212

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lysosomal acid lipase/cholesteryl ester hydrolase is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	9G7F12/7G6D7
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide.
Isotype	IgG2a
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	3988
Gene Symbol	LIPA
Uniprot ID	LICH_HUMAN
Immunogen	Purified recombinant fragment of LAL expressed in E. Coli.
Immunogen Region	
Specificity	LAL Monoclonal Antibody detects endogenous levels of LAL protein.
Immunogen Sequence	



Western blot analysis using LAL monoclonal antibody against LAL recombinant protein.