

## Anti-STK11 antibody [4H12] (STJ98220)

STJ98220

### GENERAL INFORMATION



<b>Product Type</b>	Primary antibodies
<b>Short Description</b>	Mouse monoclonal antibody anti-Serine/threonine-protein kinase STK11 is suitable for use in Western Blot, Flow Cytometry and ELISA research applications.
<b>Applications</b>	WB/FC/ELISA
<b>Host/Source</b>	Mouse
<b>Reactivity</b>	Human/Mouse/Monkey

### PRODUCT PROPERTIES

<b>Clonality</b>	Monoclonal
<b>Clone ID</b>	4H12
<b>Concentration</b>	
<b>Conjugation</b>	Unconjugated
<b>Purification</b>	Affinity purification
<b>Dilution Range</b>	WB 1:500-1:2000 FC 1:200-1:400 ELISA 1:10000
<b>Formulation</b>	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
<b>Isotype</b>	IgG1
<b>Storage Instruction</b>	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

### TARGET INFORMATION

<b>Gene ID</b>	6794
<b>Gene Symbol</b>	STK11
<b>Uniprot ID</b>	STK11_HUMAN
<b>Immunogen</b>	Purified recombinant fragment of human LKB1 expressed in E. Coli.
<b>Immunogen Region</b>	
<b>Specificity</b>	LKB1 Monoclonal Antibody detects endogenous levels of LKB1 protein.
<b>Immunogen Sequence</b>	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.  
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081