

Anti-NKX2-5 antibody [2E1] (STJ98278)

STJ98278

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Homeobox protein Nkx-2.5 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2E1
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Antibody are purified by protein G affinity chromatography.
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in 0.01M PBS pH7.4 containing 0.03% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	1482
Gene Symbol	NKX2-5
Uniprot ID	NKX25_HUMAN
Immunogen	Purified recombinant fragment of human Nkx-2.5 expressed in E. Coli.
Immunogen Region	
Specificity	Nkx-2.5 Monoclonal Antibody detects endogenous levels of Nkx-2.5 protein.
Immunogen Sequence	



This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081