

Anti-OVA antibody [3G2E1D9/3D5F4A8/3G1E1D9] (STJ98294)

STJ98294

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-OVA is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3G2E1D9/3D5F4A8/3G1E1D9
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Ovalbumin
Immunogen Region	
Specificity	OVA Monoclonal Antibody detects OVA protein.
Immunogen Sequence	



Western blot analysis using OVA monoclonal antibody against OVA protein.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081