

Anti-PAX8 antibody [2D4] (STJ98311)

STJ98311

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Paired box protein Pax-8 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	2D4
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7849
Gene Symbol	PAX8
Uniprot ID	PAX8_HUMAN
Immunogen	Purified recombinant fragment of human Pax-8 expressed in E. Coli.
Immunogen Region	
Specificity	Pax-8 Monoclonal Antibody detects endogenous levels of Pax-8 protein.
Immunogen Sequence	



Western blot analysis using Pax-8 monoclonal antibody against HeLa (1), HEK293 (2) and Raji (3) cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081