

Anti-ZBTB16 antibody [5B3] (STJ98325)

STJ98325

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Zinc finger and BTB domain-containing protein 16 is suitable for use in Western Blot, Immunofluorescence and ELISA research applications.
Applications	WB/IF/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES


Clonality	Monoclonal
Clone ID	5B3
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 IF 1:200-1:1000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7704
Gene Symbol	ZBTB16
Uniprot ID	ZBT16_HUMAN
Immunogen	Purified recombinant fragment of human PLZF expressed in E. Coli.
Immunogen Region	
Specificity	PLZF Monoclonal Antibody detects endogenous levels of PLZF protein.
Immunogen Sequence	



Immunofluorescence analysis of HeLa cells using PLZF monoclonal antibody (green). Red: Actin filaments have been labeled with Alexa Fluor-555 phalloidin.



Western blot analysis using PLZF monoclonal antibody against HeLa (1) cell lysate.