

Anti-PPARG antibody (170-270 aa) [3A4A9/1E6A1] (STJ98330)

STJ98330

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Peroxisome proliferator-activated receptor gamma (170-270 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3A4A9/1E6A1
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5468
Gene Symbol	PPARG
Uniprot ID	PPARG_HUMAN
Immunogen	Purified recombinant fragment of PPAR-Gamma (aa170-270) expressed in E. Coli.
Immunogen Region	170-270 aa
Specificity	PPAR-Gamma Monoclonal Antibody detects endogenous levels of PPAR-Gamma protein.
Immunogen Sequence	



Western blot analysis using PPAR-Gamma monoclonal antibody against truncated PPARG-His recombinant protein (1) and full-length PPARG (aa1-477) transfected CHO-K1 cell lysate (2).