

Anti-Ractopamine antibody [Rac1] (STJ98349)

STJ98349

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Ractopamine is suitable for use in ELISA research applications.
Applications	ELISA
Host/Source	Mouse
Reactivity	Species independent

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	Rac1
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	
Gene Symbol	
Uniprot ID	
Immunogen	Rac-OVA.
Immunogen Region	
Specificity	Ractopamine Monoclonal Antibody detects Ractopamine protein.
Immunogen Sequence	