

Anti-TNK1 antibody (451-560 aa) [1B5G3/3F11H6] (STJ98425)

STJ98425

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Non-receptor tyrosine-protein kinase TNK1 (451-560 aa) is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1B5G3/3F11H6
Concentration	
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:500-1:2000 ELISA 1:10000
Formulation	Liquid in PBS containing 0.03% Sodium Azide, 0.5% BSA, 50% Glycerol.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8711
Gene Symbol	TNK1
Uniprot ID	TNK1_HUMAN
Immunogen	Purified recombinant fragment of TNK1 (aa451-560) expressed in E. Coli.
Immunogen Region	451-560 aa
Specificity	TNK1 Monoclonal Antibody detects endogenous levels of TNK1 protein.
Immunogen Sequence	



Western blot analysis using TNK1 monoclonal antibody against truncated TNK1-His recombinant protein (1).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081