

Anti-IKZF1 antibody (C-Term) (STJ98508)

STJ98508

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-DNA-binding protein Ikaros (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Dog

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	10320
Gene Symbol	IKZF1
Uniprot ID	IKZF1_HUMAN
Immunogen	Purified recombinant human Ikaros (C-terminus) protein fragments expressed in E.coli.
Immunogen Region	C-Term
Specificity	Ikaros Monoclonal Antibody detects endogenous levels of Ikaros protein.
Immunogen Sequence	



Western blot analysis using Ikaros monoclonal antibody against Ramos, Jurkat cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081