

Anti-PCK2 antibody (C-Term) (STJ98530)

STJ98530

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Phosphoenolpyruvate carboxykinase, mitochondrial (C-Term) is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat/Cow/Dog/Pig

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Isotype	
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	5106
Gene Symbol	PCK2
Uniprot ID	PCKGM_HUMAN
Immunogen	Purified recombinant human PEPCK (C-terminus) protein fragments expressed in E.coli.
Immunogen Region	C-Term
Specificity	PEPCK Monoclonal Antibody detects endogenous levels of PEPCK protein.
Immunogen Sequence	



Western blot analysis using PEPCK monoclonal antibody against ZR75-1, HepG2, HeLa cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081