

Anti-SAFB antibody (STJ98547)

STJ98547

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Scaffold attachment factor B1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	6294
Gene Symbol	SAFB
Uniprot ID	SAFB1_HUMAN
Immunogen	Purified recombinant human SAF-B protein fragments expressed in E.coli.
Immunogen Region	
Specificity	SAF-B Monoclonal Antibody detects endogenous levels of SAF-B protein.
Immunogen Sequence	



Western blot analysis using SAF-B monoclonal antibody against HeLa nuclear extract, A431, MDA-MB-468, breast cancer, ZR75-1, MCF7 cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081