

Anti-WRAP53 antibody (STJ98569)

STJ98569

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Telomerase Cajal body protein 1 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	Affinity purification
Dilution Range	WB 1:1000-1:2000
Formulation	Liquid in 0.1M Tris-Glycine (pH7.4, 150 mM NaCl) , 0.2% Sodium Azide and 50% Glycerol.
Isotype	
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	55135
Gene Symbol	WRAP53
Uniprot ID	TCAB1_HUMAN
Immunogen	Purified recombinant human WDR79 protein fragments expressed in E.coli.
Immunogen Region	
Specificity	WDR79 Monoclonal Antibody detects endogenous levels of WDR79 protein.
Immunogen Sequence	



Western blot analysis using WDR79 monoclonal antibody against HeLa, U2OS cell lysate.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081