

Anti-TIMP4 antibody (174-223 aa) (STJ98630)

STJ98630

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Metalloproteinase inhibitor 4 (174-223 aa) is suitable for use in Western Blot, ELISA and Immunohistochemistry research applications.
Applications	WB/ELISA/IHC
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500-2000 IHC-P 1:50-300 ELISA 2000-20000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	7079
Gene Symbol	TIMP4
Uniprot ID	TIMP4_HUMAN
Immunogen	Synthesized peptide derived from TIMP-4 at the amino acid range 174-223
Immunogen Region	174-223 aa
Specificity	TIMP-4 Polyclonal Antibody detects endogenous levels of TIMP-4
Immunogen Sequence	



Western blot analysis of 3T3 cells using TIMP-4 Polyclonal Antibody diluted at 1:800. Secondary antibody was diluted at 1:20000