

Anti-BTN3A1/2/3 antibody (41-90 aa) (STJ98764)

STJ98764

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Butyrophilin subfamily 3 member A1 and Butyrophilin subfamily 3 member A2 and Butyrophilin subfamily 3 member A3 (41-90 aa) is suitable for use in Immunohistochemistry, Immunofluorescence and ELISA research application
Applications	IHC/IF/ELISA
Host/Source	Rabbit
Reactivity	Human/Rat/Mouse

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution	IHC-P 1:50-200
Range	ELISA 1:10000-20000 IF 1:50-200
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	11119 11118 10384 BTN3A1 BTN3A2 BT3A1_HUMAN BT3A2_HUMAN BT3A3_HUMAN
Immunogen	Synthetic peptide from the human protein at the amino acid range 41-90
Immunogen Region	41-90 aa
Specificity	The antibody detects endogenous BTN3A1/2/3
Immunogen Sequence	



Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

Immunohistochemical analysis of paraffin-embedded human-kidney, antibody was diluted at 1:200

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081