

Anti-Di-Methyl-Histone H3-K80 antibody (STJ98849)

STJ98849

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Rabbit polyclonal antibody anti-Di-Methyl-Histone H3.1 and Histone H3.2 and Histone H3.3 and Histone H3.3C-K80 is suitable for use in Western Blot and ELISA research applications.
Applications	WB/ELISA
Host/Source	Rabbit
Reactivity	Human/Mouse/Rat

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	440093 3020/3021
Gene Symbol	H3-5 H3-3A.H3-3B H3C_HUMAN H33_HUMAN H32_HUMAN <
Immunogen	Synthesized Tri-Methyl peptide derived from the human Histone H3. at the amino acid range K80
Immunogen Region	
Specificity	This antibody detects endogenous levels of Histone H3 at Human:K80Mouse:K80Rat:K80, It doesn't react with total protein.
Immunogen Sequence	



Western blot analysis of HeLa cells using primary antibody diluted at 1:1000 (4A°C overnight). Secondary antibody: IgG-coated Anti-rabbit IgG IRDye 800 (diluted at 1:5000, 25A°C, 1 hour)

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081