

Anti-HNRNPU antibody [1B11-H2-B7] (STJ99127)

STJ99127

GENERAL INFORMATION

Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Heterogeneous nuclear ribonucleoprotein U is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	1B11-H2-B7
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	3192
Gene Symbol	HNRNPU
Uniprot ID	HNRPU_HUMAN
Immunogen	Purified recombinant human hnRNP U protein fragments expressed in E.coli.
Immunogen Region	
Specificity	This antibody detects endogenous levels of hnRNP U and does not cross-react with related proteins.
Immunogen Sequence	



Western blot detection of hnRNP U in K562, Jurkat and HeLa cell lysates and using hnRNP U mouse mAb (1:1000 diluted). Predicted band size: 110 kDa. Observed band size: 110 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081