

Anti-ANXA6 antibody [8A7-F8-G12] (STJ99178)

STJ99178

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Annexin A6 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8A7-F8-G12
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	309
Gene Symbol	ANXA6
Uniprot ID	ANXA6_HUMAN
Immunogen	Purified recombinant human Annexin VI protein fragments expressed in E.coli.
Immunogen Region	
Specificity	This antibody detects endogenous levels of Annexin VI and does not cross-react with related proteins.
Immunogen Sequence	



Western blot detection of Annexin VI in T47D, MB4, CEM, MOLT-4 and Jurkat cell lysates and using Annexin VI mouse mAb (1:1000 diluted). Predicted band size: 76 kDa. Observed band size: 76 kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081