

Anti-EIF4EBP3 antibody [3B6-D3-H7] (STJ99277)

STJ99277

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Eukaryotic translation initiation factor 4E-binding protein 3 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Transfected

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	3B6-D3-H7
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.
Instruction	

TARGET INFORMATION

Gene ID	8637
Gene Symbol	EIF4EBP3
Uniprot ID	4EBP3_HUMAN
Immunogen	Recombinant protein of human EIF4EBP3.
Immunogen	
Region	
Specificity	Transfected Only.
Immunogen	
Sequence	



Western blot detection of EIF4EBP3 in CHO-K1 cell lysate (B) and CHO-K1 transfected by EGFP-EIF4EBP3 fragment (A) cell lysate using EIF4EBP3 mouse mAb (1:1000 diluted).

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081