

Anti-LMNB2 antibody [8F6-E8-F12] (STJ99315)

STJ99315

GENERAL INFORMATION


Product Type	Primary antibodies
Short Description	Mouse monoclonal antibody anti-Lamin-B2 is suitable for use in Western Blot research applications.
Applications	WB
Host/Source	Mouse
Reactivity	Human/Mouse

PRODUCT PROPERTIES

Clonality	Monoclonal
Clone ID	8F6-E8-F12
Concentration	1 mg/mL
Conjugation	Unconjugated
Purification	The antibody was affinity-purified from mouse ascites by affinity-chromatography using epitope-specific immunogen.
Dilution Range	WB 1:500
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG1
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID	84823
Gene Symbol	LMNB2
Uniprot ID	LMNB2_HUMAN
Immunogen	Recombinant human Lamin B2 protein.
Immunogen Region	
Specificity	This antibody detects endogenous levels of Lamin B2 and does not cross-react with related proteins.
Immunogen Sequence	



Western blot detection of Lamin B2 in HCT116, Mouse Embryo and Jurkat cell lysates using Lamin B2 mouse mAb (dilution 1:500). Predicted band size: 68kDa. Observed band size: 68kDa.

This product is suitable for in-vitro studies under the RESEARCH USE ONLY [RUO] licence. This product must not be used as for diagnostic or other medical purposes.
St John's Laboratory Ltd, Knowledge Dock Business Centre, University Way, London, E16 2RD | Tel: 0208 223 3081