

Goat Anti-Rabbit IgG Heavy & Light Chain antibody {Gold} (STJ99356)

STJ99356

GENERAL INFORMATION

Product Type	Secondary antibodies
Short Description	Gold-conjugated goat polyclonal anti-Rabbit IgG Heavy & Light Chain secondary antibody. For use in most research applications.
Applications	WB/IHC/IF/ELISA
Host/Source	Goat
Reactivity	Rabbit

PRODUCT PROPERTIES

Clonality	Polyclonal
Clone ID	
Concentration	2 mg/mL
Conjugation	Gold
Purification	The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen.
Dilution Range	WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000
Formulation	Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide.
Isotype	IgG
Storage Instruction	Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles.

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen
Sequence