

Goat Anti-Tree Shrew IgG Heavy & Light Chain antibody {Gold} (STJ99388)

STJ99388

GENERAL INFORMATION

| | |
|--------------------------|--|
| Product Type | Secondary antibodies |
| Short Description | Gold-conjugated goat polyclonal anti-Tree Shrew IgG Heavy & Light Chain secondary antibody. For use in most research applications. |
| Applications | WB/IHC/IF/ELISA |
| Host/Source | Goat |
| Reactivity | Tree Shrew |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Polyclonal |
| Clone ID | |
| Concentration | 2 mg/mL |
| Conjugation | Gold |
| Purification | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range | WB 1:1000-10000 |
| | ELISA 1:10000-20000 |
| | IHC 1:50-100 |
| | IF 1:50-1000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

| |
|-------------------------|
| Gene ID |
| Gene Symbol |
| Uniprot ID |
| Immunogen |
| Immunogen Region |
| Specificity |
| Immunogen |
| Sequence |