

Mouse Anti-Human IgG3 antibody {Gold} (STJ99394)

STJ99394

GENERAL INFORMATION

| | |
|--------------------------|---|
| Product Type | Secondary antibodies |
| Short Description | Gold-conjugated mouse monoclonal anti-Human IgG3 secondary antibody. For use in most research applications. |
| Applications | WB/IHC/IF/ELISA |
| Host/Source | Mouse |
| Reactivity | Human |

PRODUCT PROPERTIES

| | |
|----------------------------|---|
| Clonality | Monoclonal |
| Clone ID | |
| Concentration | 2 mg/mL |
| Conjugation | Gold |
| Purification | The antibody was affinity-purified by affinity-chromatography using antigen-specific immunogen. |
| Dilution Range | WB 1:1000-10000 ELISA 1:10000-20000 IHC 1:50-100 IF 1:50-1000 |
| Formulation | Liquid in PBS containing 50% Glycerol, 0.5% BSA and 0.02% Sodium Azide. |
| Isotype | IgG3 |
| Storage Instruction | Store at -20°C for up to 1 year from the date of receipt, and avoid repeat freeze-thaw cycles. |

TARGET INFORMATION

Gene ID
Gene Symbol
Uniprot ID
Immunogen
Immunogen Region
Specificity
Immunogen Sequence